

Electric Patient Lift

ML-10



Assisting device for the disabled person

Man&tel
주/맨엔텔

Man&tel Co., Ltd.
258 Saneop-ro, Gumi-si, Gyeongsangbuk-do 39232, Korea
Tel. +82-54-444-2377
www.manntel.com master@manntel.com

Electric patient lift

ML-10

- Item : Electric patient lift
- Patient to use this device: Patient with discomfort in the behavior and movement, disable person

Man&tel
주/맨엔텔

Assisting device for the disabled person Item of health insurance benefit

Electric patient lift giving care partner convenience and stability to the patient when patient move to wheel chair, bed and bathroom from the floor
Convenient device for the physically challenged, hemiplegic patient with cervical spine disorders, with thoracic spine disorders, and with mobility difficulties and for the care partners

Medical equipment certification

- KFDA Product name : Lift, electrically-powered
- KFDA Free sales Certification No. 16-646
- KFDA Manufacturing Permission No. 3399
- KGMP Certificate
- ISO13485 Certificate
- Korea assisting device items management No. MA0012000011

Features

■ Basic

- Able to move patient to bed and wheel chair from floor
- Full charge able to use 75 times
- When using charging system, not needed AC power
- Very strong steeliness(Safe weight loading 150kg)
- Height adjustable using remote control unit
- Adjustable range of height to move patient(0~770mm)
- Dual purpose lifting sling to use for cervical vertebral
- Folded type for easy storage and move

■ Option

- Measurement of weight

Product specification

- Maximum patient weight: 150Kg
- Product size : 650(W) X 1,200(L) X 1,900(H)mm ±5%
- Weight: 55Kg
- Control box Adaptor INPUT : AC 220~240V, 50/60Hz, 58.8W
- Control box Adaptor OUTPUT : DC 29V, 2A
- Control box Battery : lithium polymer DC 24V, 4.4[AH]
- Weight machine(option) Battery : 1.5V(AAA) X 2

Composition

- Control box, Cylinder, Battery, Remote control unit, Adaptor, Patient protection lifting sling, Weight machine(option)



※ This device is a medical equipment. Prior to use, please read carefully precautions to use and the way of use.

Published: Oct, 2019

www.manntel.com
Tel. 054-444-2377, Email. master@manntel.com